

Factor IX In Vitro Transcribed mRNA-LNP

Catalog Number:SG-MRNA-LNP-1871

DESCRIPTION	
Product Name	Factor IX In Vitro Transcribed mRNA-LNP
Gene Name	Factor IX
Source	The ORF of Factor IX was cloned in our IVT vector and mRNA was prepared through in vitro transcription and purification. The purified mRNA was further encapsulated with LNP(DSPC:Cholesterol:DMG-PEG:SM102).
Alternative names	Factor IX
SPECIFICATIONS	
Cap	m7GpppN
5'-UTR	5' -untranslated region derived from human alpha-globin RNA with an optimized Kozak sequence
ORF	Factor IX
3'-UTR	3' UTR comprising two sequence elements derived from the aminoterminal enhancer of split (AES) mRNA and the mitochondrial encoded 12S ribosomal RNA
Poly(A) Tail	A 110-nucleotide poly(A)-tail consisting of a stretch of 30 adenosine residues, followed by a 10-nucleotide linker sequence and another 70 adenosine residues.
Modifications	N1-methyl-pseudouridine
Neutral Lipid	1,2-distearoyl-sn-glycero-3-phosphocholine (DSPC)
Cholesterol	Cholesterol
Ionizable Lipid	1,2-dimyristoyl-rac-glycero-3-methoxypolyethylene glycol-2000 (PEG2000-DMG)
PEG-lipid	Heptadecan-9-yl 8-((2-hydroxyethyl)(8-(nonyloxy)-8-oxooctyl)amino)octanoate)(SM-102)
Storage	-20 °C
Buffer	PBS, pH7.4
Cryoprotectant	Trehalose
BACKGROUND	
Gene Accession	
Gene Alias	Factor IX
Background	Factor IX (or Christmas factor) is one of the serine proteases of the coagulation system; it belongs to peptidase family S1. Deficiency of this protein causes hemophilia B.



Seattle Genova

Tel: +1 (425) 247-3088 Fax: +1 (425) 650-9990

Email: info@seattle-genova.com Web: www.seattle-genova.com

Address: 18110 SE 34TH ST STE 455, Vancouver, WA 98683
