

Human deoxyribonuclease I In Vitro Transcribed mRNA-LNP

Catalog Number:SG-MRNA-LNP-1918

DESCRIPTION	
Product Name	Human deoxyribonuclease I In Vitro Transcribed mRNA-LNP
Gene Name	DNase I
Source	The ORF of Human deoxyribonuclease I was cloned in our IVT vector and mRNA was prepared through in vitro transcription and purification. The purified mRNA was further encapsulated with LNP(DSPC:Cholesterol:DMG-PEG:SM102).
Alternative names	Human deoxyribonuclease I
SPECIFICATIONS	
Cap	m7GpppN
5'-UTR	5' -untranslated region derived from human alpha-globin RNA with an optimized Kozak sequence
ORF	Human deoxyribonuclease I
3'-UTR	3' UTR comprising two sequence elements derived from the aminoterminal enhancer of split (AES) mRNA and the mitochondrial encoded 12S ribosomal RNA
Poly(A) Tail	A 110-nucleotide poly(A)-tail consisting of a stretch of 30 adenosine residues, followed by a 10-nucleotide linker sequence and another 70 adenosine residues.
Modifications	N1-methyl-pseudouridine
Neutral Lipid	1,2-distearoyl-sn-glycero-3-phosphocholine (DSPC)
Cholesterol	Cholesterol
Ionizable Lipid	1,2-dimyristoyl-rac-glycero-3-methoxypolyethylene glycol-2000 (PEG2000-DMG)
PEG-lipid	Heptadecan-9-yl 8-((2-hydroxyethyl)(8-(nonyloxy)-8-oxooctyl)amino)octanoate)(SM-102)
Storage	-20 °C
Buffer	PBS, pH7.4
Cryoprotectant	Trehalose
BACKGROUND	
Gene Accession	
Gene Alias	Human deoxyribonuclease I
Background	Human deoxyribonuclease I (DNase I) is an endonuclease that catalyzes the hydrolysis of extracellular DNA and is just one of the numerous types of nucleases found in nature.



Seattle Genova

Tel: +1 (425) 247-3088 Fax: +1 (425) 650-9990

Email: info@seattle-genova.com Web: www.seattle-genova.com

Address: 18110 SE 34TH ST STE 455, Vancouver, WA 98683
