

Enhanced yellow fluorescent protein(EYFP) In Vitro Transcribed mRNA-LNP

Catalog Number:MRNA-TG-016

DESCRIPTION	
Product Name	Enhanced yellow fluorescent protein(EYFP) In Vitro Transcribed mRNA-LNP
Gene Name	Enhanced yellow fluorescent protein(EYFP)
Source	In vitro transcribed mRNA encapsulated with LNP
Alternative names	
SPECIFICATIONS	
Cap	Cap 1
5'-UTR	5' -untranslated region derived from human alpha-globin RNA with an optimized Kozak sequence
ORF	Enhanced yellow fluorescent protein(EYFP)
3'-UTR	3' UTR comprising two sequence elements derived from the aminoterminal enhancer of split (AES) mRNA and the mitochondrial encoded 12S ribosomal RNA
Poly(A) Tail	A 110-nucleotide poly(A)-tail consisting of a stretch of 30 adenosine residues, followed by a 10-nucleotide linker sequence and another 70 adenosine residues.
Modifications	N1-methyl-pseudouridine
Neutral Lipid	1,2-distearoyl-sn-glycero-3-phosphocholine (DSPC)
Cholesterol	Cholesterol
Ionizable Lipid	1,2-dimyristoyl-rac-glycero-3-methoxypolyethylene glycol-2000 (PEG2000-DMG)
PEG-lipid	Heptadecan-9-yl 8-((2-hydroxyethyl)(8-(nonyloxy)-8-oxooctyl)amino)octanoate)(SM-102)
Storage	-20 °C
Buffer	PBS, pH7.4
Cryoprotectant	Trehalose
BACKGROUND	
Gene Accession	
Gene Alias	
Background	EYFP is a basic (constitutively fluorescent) green/yellow fluorescent protein published in 1996, derived from <i>Aequorea victoria</i> . It is reported to be a very rapidly-maturing weak dimer with high acid sensitivity.



Seattle Genova

Tel: +1 (425) 247-3088 Fax: +1 (425) 650-9990

Email: info@seattle-genova.com Web: www.seattle-genova.com

Address: 18110 SE 34TH ST STE 455, Vancouver, WA 98683
