

Aed a 3 Circular RNA

Catalog Number:ALG-ORNA-0026

DESCRIPTION	
Product Name	Aed a 3 Circular RNA
Gene Name	Aed a 3
Source	In vitro transcribed mRNA was further circularized to make this product as a circular RNA.
Alternative names	Uncharacterized 30 kDa salivary protein
SPECIFICATIONS	
Cap	
5'-UTR	
ORF	Aed a 3
3'-UTR	
Poly(A) Tail	
Modifications	
Neutral Lipid	1,2-distearoyl-sn-glycero-3-phosphocholine (DSPC)
Cholesterol	Cholesterol
Lonizable Lipid	1,2-dimyristoyl-rac-glycero-3-methoxypolyethylene glycol-2000 (PEG2000-DMG)
PEG-lipid	Heptadecan-9-yl 8-((2-hydroxyethyl)(8-(nonyloxy)-8-oxooctyl)amino)octanoate)(SM-102)
Storage	-80 °C
Buffer	PBS, pH7.30
Cryoprotectant	Trehalose
BACKGROUND	
Gene Accession	
Gene Alias	Uncharacterized 30 kDa salivary protein
Background	Of 28 mosquito-bite test-positive subjects, 9 (32%) had a positive skin reaction to rAed a 3.