

Aed al 2 Circular RNA

Catalog Number:ALG-ORNA-0034

DESCRIPTION	
Product Name	Aed al 2 Circular RNA
Gene Name	Aed al 2
Source	In vitro transcribed mRNA was further circularized to make this product as a circular RNA.
Alternative names	D7 like salivary odorant binding protein
SPECIFICATIONS	
Cap	
5'-UTR	
ORF	Aed al 2
3'-UTR	
Poly(A) Tail	
Modifications	
Neutral Lipid	1,2-distearoyl-sn-glycero-3-phosphocholine (DSPC)
Cholesterol	Cholesterol
Lonizable Lipid	1,2-dimyristoyl-rac-glycero-3-methoxypolyethylene glycol-2000 (PEG2000-DMG)
PEG-lipid	Heptadecan-9-yl 8-((2-hydroxyethyl)(8-(nonyloxy)-8-oxooctyl)amino)octanoate)(SM-102)
Storage	-80 °C
Buffer	PBS, pH7.38
Cryoprotectant	Trehalose
BACKGROUND	
Gene Accession	
Gene Alias	D7 like salivary odorant binding protein
Background	Allergenicity reference:PubMed:32099589 Route of allergen exposure:Injection; 10 of 15 subjects with immediate or delayed skin reactivity to bites by Asian Tiger mosquito were positive in IgE binding to this protein by western blot.