

Alt a 5 Circular RNA

Catalog Number:ALG-ORNA-0047

DESCRIPTION	
Product Name	Alt a 5 Circular RNA
Gene Name	Alt a 5
Source	In vitro transcribed mRNA was further circularized to make this product as a circular RNA.
Alternative names	Ribosomal protein P2
SPECIFICATIONS	
Cap	
5'-UTR	
ORF	Alt a 5
3'-UTR	
Poly(A) Tail	
Modifications	
Neutral Lipid	1,2-distearoyl-sn-glycero-3-phosphocholine (DSPC)
Cholesterol	Cholesterol
Lonizable Lipid	1,2-dimyristoyl-rac-glycero-3-methoxypolyethylene glycol-2000 (PEG2000-DMG)
PEG-lipid	Heptadecan-9-yl 8-((2-hydroxyethyl)(8-(nonyloxy)-8-oxooctyl)amino)octanoate)(SM-102)
Storage	-80 °C
Buffer	PBS, pH7.51
Cryoprotectant	Trehalose
BACKGROUND	
Gene Accession	
Gene Alias	Ribosomal protein P2
Background	Obsolete name:Alt a 6; Allergenicity:- 6 out of 43 sera (14%) from A.alternaria-sensitive patients showed IgE binding to rAlt a 6 on immunoblot. - Of 98 A.alternata-sensitized patients (positive immunoblot of A.alternata extract), 8% showed IgE binding to Alt a 6 on immunoblot: medline 95206305 and 95337840