

Alt a 6 Circular RNA

Catalog Number:ALG-ORNA-0048

DESCRIPTION	
Product Name	Alt a 6 Circular RNA
Gene Name	Alt a 6
Source	In vitro transcribed mRNA was further circularized to make this product as a circular RNA.
Alternative names	Enolase
SPECIFICATIONS	
Cap	
5'-UTR	
ORF	Alt a 6
3'-UTR	
Poly(A) Tail	
Modifications	
Neutral Lipid	1,2-distearoyl-sn-glycero-3-phosphocholine (DSPC)
Cholesterol	Cholesterol
Lonizable Lipid	1,2-dimyristoyl-rac-glycero-3-methoxypolyethylene glycol-2000 (PEG2000-DMG)
PEG-lipid	Heptadecan-9-yl 8-((2-hydroxyethyl)(8-(nonyloxy)-8-oxooctyl)amino)octanoate)(SM-102)
Storage	-80 °C
Buffer	PBS, pH7.52
Cryoprotectant	Trehalose
BACKGROUND	
Gene Accession	
Gene Alias	Enolase
Background	Obsolete name:Alt a 5, 11; Allergenicity:- 2 out of 7 A.alternaria-sensitive patients gave positive skin reactions in response to Alt a 11. - Of 23 A.alternaria-sensitized patients tested, 22% showed serum IgE binding to purified Alt a 11 on immunoblot: medline 20534817.