

Amb a 11 Circular RNA

Catalog Number:ALG-ORNA-0055

DESCRIPTION	
Product Name	Amb a 11 Circular RNA
Gene Name	Amb a 11
Source	In vitro transcribed mRNA was further circularized to make this product as a circular RNA.
Alternative names	Cysteine protease
SPECIFICATIONS	
Cap	
5'-UTR	
ORF	Amb a 11
3'-UTR	
Poly(A) Tail	
Modifications	
Neutral Lipid	1,2-distearoyl-sn-glycero-3-phosphocholine (DSPC)
Cholesterol	Cholesterol
Ionizable Lipid	1,2-dimyristoyl-rac-glycero-3-methoxypolyethylene glycol-2000 (PEG2000-DMG)
PEG-lipid	Heptadecan-9-yl 8-((2-hydroxyethyl)(8-(nonyloxy)-8-oxooctyl)amino)octanoate)(SM-102)
Storage	-80 °C
Buffer	PBS, pH7.59
Cryoprotectant	Trehalose
BACKGROUND	
Gene Accession	
Gene Alias	Cysteine protease
Background	Allergenicity reference:PubMed:25865353 Route of allergen exposure:Airway; Sera from 38 of 70 short ragweed pollen allergic individuals contained IgE that reacted with nAmb a 11 on 2D immunoblots.