

## Ani s 13 Circular RNA

Catalog Number:ALG-ORNA-0078

DESCRIPTION	
Product Name	Ani s 13 Circular RNA
Gene Name	Ani s 13
Source	In vitro transcribed mRNA was further circularized to make this product as a circular RNA.
Alternative names	Hemoglobin
SPECIFICATIONS	
Cap	
5'-UTR	
ORF	Ani s 13
3'-UTR	
Poly(A) Tail	
Modifications	
Neutral Lipid	1,2-distearoyl-sn-glycero-3-phosphocholine (DSPC)
Cholesterol	Cholesterol
Ionizable Lipid	1,2-dimyristoyl-rac-glycero-3-methoxypolyethylene glycol-2000 (PEG2000-DMG)
PEG-lipid	Heptadecan-9-yl 8-((2-hydroxyethyl)(8-(nonyloxy)-8-oxooctyl)amino)octanoate)(SM-102)
Storage	-80 °C
Buffer	PBS, pH7.82
Cryoprotectant	Trehalose
BACKGROUND	
Gene Accession	
Gene Alias	Hemoglobin
Background	Allergenicity reference:PubMed:25683373 Route of allergen exposure:Food; Twenty eight of 44 case-history positive subjects diagnosed with Gastro-allergic Anisakiasis or Anisakis sensitization associated chronic urticaria, had demonstrable IgE binding to the protein by hemoglobin-specific capture ELISA.