

Vig r 6 Circular RNA

Catalog Number:ALG-ORNA-1001

DESCRIPTION	
Product Name	Vig r 6 Circular RNA
Gene Name	Vig r 6
Source	In vitro transcribed mRNA was further circularized to make this product as a circular RNA.
Alternative names	Cytokinin-specific binding protein (CSBP), Bet v 1 family member
SPECIFICATIONS	
Cap	
5'-UTR	
ORF	Vig r 6
3'-UTR	
Poly(A) Tail	
Modifications	
Neutral Lipid	1,2-distearoyl-sn-glycero-3-phosphocholine (DSPC)
Cholesterol	Cholesterol
Ionizable Lipid	1,2-dimyristoyl-rac-glycero-3-methoxypolyethylene glycol-2000 (PEG2000-DMG)
PEG-lipid	Heptadecan-9-yl 8-((2-hydroxyethyl)(8-(nonyloxy)-8-oxooctyl)amino)octanoate)(SM-102)
Storage	-80 °C
Buffer	PBS, pH7.1005
Cryoprotectant	Trehalose
BACKGROUND	
Gene Accession	
Gene Alias	Cytokinin-specific binding protein (CSBP), Bet v 1 family member
Background	Allergenicity reference:PubMed:23996905 Route of allergen exposure:Food; Of 19 patients who reported allergic reactions or had positive prick-to-prick tests to mung bean sprouts, 63% showed IgE binding to Vig r 6 by ELISA.